

## BOTTLES - "for water sampling" - "PET"



➤ Manufactured from transparent PET with thick walls and crush-resistant corners.

➤ Ideal for taking, transporting and storing of samples for microbiological analysis.

➤ They can be used for all analysis where a sterile sample is required.

➤ Crystall clear bottle body enables the user to visibly check the sample integrity.

➤ All other technical and physical specs are similar as P.P bottles.

bulk packed bottles	bottle volume	bottle dimensions	clear bottles W/O sodium thiosulfate	clear bottles with sodium thiosulfate	packing specifications
	125 ml	56 x 56 x 90 mm	061.61.125	061.63.125	180 pcs in box
250 ml	60 x 60 x 116 mm	061.61.250	061.63.250	108 pcs in box	
500 ml	70 x 70 x 145 mm	061.61.500	061.63.500	84 pcs in box	
1000 ml	90 x 90 x 174 mm	061.61.901	061.63.901	30 pcs in box	

single packed bottles	bottle volume	bottle dimensions	clear bottles W/O sodium thiosulfate	clear bottles with sodium thiosulfate	packing specifications
	125 ml	56 x 56 x 90 mm	061.65.125	061.67.125	130 pcs in box
250 ml	60 x 60 x 116 mm	061.65.250	061.67.250	72 pcs in box	
500 ml	70 x 70 x 145 mm	061.65.500	061.67.500	42 pcs in box	
1000 ml	90 x 90 x 174 mm	061.65.901	061.67.901	22 pcs in box	



## CONTAINER CUP

➤ Manufactured from clear polypropylene with 40 ml effective volume.

➤ The press-on cap ensures hermetical and leak-proof closure. Scale on side wall for easy volume identification.

➤ Frosted marking area on side wall enables easy marking and labelling.

JARS  
"glass" - "wide neck"

➤ Manufactured from neutral glass with thick walls and crush-resistant corners.

➤ Ideal for safe storage and transportation of liquids and powders.

➤ Enamelled steel cap inner design ensures hermetical and leak-proof closure.

catalogue number	product description	pack quantity
087.02.001	non-sterile cups	1000 pieces
087.02.002	aseptic sterile cups	500 pieces

catalogue number	volume	jar diameter	pack quantity
061.29.040	40 ml	43 mm	10 pieces
061.29.105	105 ml	49 mm	10 pieces
061.29.370	370 ml	72 mm	10 pieces
061.29.720	720 ml	90 mm	10 pieces
061.29.901	1000 ml	95 mm	10 pieces