



BOTTLES - "spray"

- Manufactured from transparent PET with a spray head.
- All media contacting parts perfectly resist to most solvents and chemicals.
- Liquids can be finely vaporised with the pull of the hand pump trigger.
- Feature with a safety lock at the trigger which prevents accidental use.

| catalogue number | volume | pack quantity |
|------------------|--------|---------------|
| 062.11.250 | 250 ml | 1 piece |
| 062.11.500 | 500 ml | 1 piece |

Extremely fine and homogenous vaporisation is guaranteed with the pull of the trigger.



WASH BOTTLES - "integral"

- Manufactured from highly translucent polyethylene. Ergonomic shape offers easy grip and easy squeeze during use.
- Integral moulded delivery tube delivers the liquid in the bottle until the last drop.
- Delivery tube guarantees a fine and straight flow with an application of a slight pressure on the body.

| catalogue number | volume | pack quantity |
|------------------|--------|---------------|
| 062.08.250 | 250 ml | 1 piece |
| 062.08.500 | 500 ml | 1 piece |



Dispensing tip provides extra fine stream. The tip of the delivery tube can be cut in order to increase the liquid flow rate.



WASH BOTTLES - "narrow neck"

- Manufactured from highly translucent and semi-rigid polyethylene with thick walls and crush-resistant corners.
- Round bottle form allows un-restricted liquid flow-out without leaving any residue at the bottle neck.
- Flexible delivery tube guarantees a fine and straight flow with an application of a slight pressure on the body.
- "Goose neck" design of the delivery tube guarantees easy and precise reaching to the flask corners.
- Dispensing tip provides extra fine stream.
- Delivery tube tip can be cut to increase the flow rate.

Colored body and caps enables fast, and easy identification of bottle content.



| clear bottles | volume | coloured bottles | volume | bottle color | pack quantity |
|---------------|---------|------------------|--------|--------------|---------------|
| 062.05.100 | 100 ml | 062.05.05R | 500 ml | red | 1 piece |
| 062.05.250 | 250 ml | 062.05.05Y | 500 ml | yellow | 1 piece |
| 062.05.500 | 500 ml | 062.05.05G | 500 ml | green | 1 piece |
| 062.05.901 | 1000 ml | 062.05.05B | 500 ml | blue | 1 piece |

