

TUBE RACKS - "polypropylene" - "foldable"



➤ Manufactured from modified polypropylene which enables the racks to be used in water baths without floating.

➤ Guarantees temperature resistance between -20°C to +121°C.

➤ 4 different color options offer visual sorting and organising of different samples.

➤ Raised tube sockets guarantee firm hold of the tubes at their exact positions.

➤ The alpha numeric grid reference on the top tier of the rack ensures an easy tube location and sample positioning.

➤ Racks are supplied assembled and folded. Can be ready for use in few seconds with a simple action.

➤ These racks can be folded after use for reducing storage area.

| for tubes | well quantity | tube array | rack L x W x H | white racks catalogue no | blue racks catalogue no | red racks catalogue no | orange racks catalogue no | pack quantity |
|-----------|---------------|------------|----------------|--------------------------|-------------------------|------------------------|---------------------------|---------------|
| 13 mm | 50 | 5 x 10 | 140x275x70 mm | 079.03.001W | 079.03.001B | 079.03.001R | 079.03.001O | 1 piece |
| 16 mm | 50 | 5 x 10 | 140x275x70 mm | 079.03.002W | 079.03.002B | 079.03.002R | 079.03.002O | 1 piece |
| 18 mm | 50 | 5 x 10 | 140x275x70 mm | 079.03.003W | 079.03.003B | 079.03.003R | 079.03.003O | 1 piece |
| 20 mm | 45 | 5 x 9 | 140x275x70 mm | 079.03.004W | 079.03.004B | 079.03.004R | 079.03.004O | 1 piece |
| 30 mm | 26 | 4 x 6/7 | 140x275x70 mm | 079.03.005W | 079.03.005B | 079.03.005R | 079.03.005O | 1 piece |

