

SPATULAS - "polypropylene" - "macro spoon"



Flexible body form offers best sampling and collection options even for wet pastes.



➤ Manufactured from high quality flexible polypropylene.

➤ Most commonly used spatula. Featured with a macro spoon form for small, medium and large size particles at one end and a straight/flat spatula form at the other end.

➤ Robust and seamless construction guarantees long life cycle. Smooth polished surfaces provide fast and easy cleaning.

➤ Also offers practical solutions for mixing and stirring solutions manually by hand.

➤ Flat end has 35 degree sharp edges for easy removal of the residues on flat surfaces.

➤ Flexible body form offers best sampling and collecting.


➤ **Alternative 2 spatulas:** ➤ **Non-sterile spatulas**
➤ **Gamma sterile spatulas**

➤ Gamma sterilised sterile spatulas carry such informations as lot number, production and expiry dates etc on the packing.

total length	non-sterile spatulas	gamma sterile spatulas	pack quantity
150 mm	047.32.150	047.33.150	1 piece
180 mm	047.32.180	047.33.180	1 piece
210 mm	047.32.210	047.33.210	1 piece



Featured with a macro spoon form for small, medium and large size particles



Flat end has 35 degree sharp edges for easy removal of the residues on flat surfaces.

SPATULAS - "sterile" - "with lid" - "polystyrene"

➤ Manufactured from high quality polystyrene which resists to most acids, bases and many commonly used solvents. Guarantees temperature resistance between -10°C to +70°C.

➤ Gamma sterile spatulas are manufactured as single use. They eliminate the need of separate sterile bags or containers in sample collection.

➤ Perfect shape and size enables efficient sampling of various powder, granulate and pasty form samples.

➤ Sharp blade penetrates directly into paper or plastic bags for optimum and quick sampling process.

➤ Easy to grip and long handle increases the distance between the hand and the sample to reduce cross contamination risks.


➤ Approximate volume values are printed on scoop handles.

➤ Sterile spatulas are featured with a transparent lid which protects the samples from dust, moisture and other effects.

➤ The transparent lid should be closed immediately after the sampling to secure a sterile and contamination-free sample.

➤ Lid is featured with a frosted write on area for practical use.

catalogue number	spatula (Ø) W	spatula L1	handle L2	pack quantity
037.36.001	13 mm	150 mm	110 mm	1 piece



Sharp blade penetrates directly into paper or plastic bags for optimum and quick sampling process.

Spatulas are featured with a lid which protects the samples from dust, moisture and other effects.



CANCELLED