


distillation & filtration & separation & sample preparation

DISTILLING LINKS - "with condenser" - "Claisen head"



► Manufactured from heat resistant, durable borosilicate glass with a liebig condenser with two NS 14/23 socket is designed to be used with distillation thermometers.

► Excellently grinded joints prevent airleaks and liquid leakage.

| catalogue number | cone NS size | socket NS size | condenser length | pack quantity |
|------------------|--------------|----------------|------------------|---------------|
| 065.57.001 | NS 14/23 | NS 14/23 | 160 mm | 1 piece |
| 065.57.002 | NS 14/23 | NS 14/23 | 250 mm | 1 piece |
| 065.57.003 | NS 29/32 | NS 14/23 | 250 mm | 1 piece |
| 065.57.004 | NS 29/32 | NS 14/23 | 400 mm | 1 piece |

DISTILLING LINKS - "with condenser & vacuum adapter" - "Claisen head"



► Manufactured from heat resistant, durable borosilicate glass with a liebig condenser with two NS 14/23 socket is designed to be used with distillation thermometers. Featured with a vacuum connecting adapter.

► Excellently grinded joints prevent airleaks and liquid leakage.

| catalogue number | cone NS size | socket NS size | condenser length | pack quantity |
|------------------|--------------|----------------|------------------|---------------|
| 065.58.001 | NS 14/23 | NS 14/23 | 160 mm | 1 piece |
| 065.58.002 | NS 14/23 | NS 14/23 | 250 mm | 1 piece |
| 065.58.003 | NS 29/32 | NS 14/23 | 250 mm | 1 piece |
| 065.58.004 | NS 29/32 | NS 14/23 | 400 mm | 1 piece |

ADAPTERS - "PTFE" - "flexible"



► Manufactured from chemically inert PTFE with excellent chemical resistance, high mechanical strength and high temperature resistance between -200°C to +260°C.

► Unique flexible body form guarantees 90° angular connections of NS joints and provides corrections of misalignment in glassware assemblies.

► Extra smooth interior finish guarantees high vacuum sealing performance without a need of any lubricant.

| catalogue number | socket NS size | cone NS size | pack quantity |
|------------------|----------------|--------------|---------------|
| 065.36.014 | NS 14/23 | NS 14/23 | 1 piece |
| 065.36.019 | NS 19/26 | NS 19/26 | 1 piece |
| 065.36.024 | NS 24/29 | NS 24/29 | 1 piece |
| 065.36.029 | NS 29/32 | NS 29/32 | 1 piece |

GLASS BEADS



► Manufactured from borosilicate glass.

► They offer practical solutions specially for distillation applications with their different diameters.

| catalogue number | diameter (mm) | pack quantity |
|------------------|---------------|---------------|
| 030.60.003 | 3,0 - 3,5 | 1 kg |
| 030.60.004 | 3,5 - 4,0 | 1 kg |
| 030.60.005 | 4,0 - 5,0 | 1 kg |
| 030.60.006 | 5,0 - 6,0 | 1 kg |